

USER MANUAL

MODEL CPC1S

Created by “Checkpoint R&D Ltd.”



Read carefully before use!

In this manual, you will learn how to use the Multiple Physiological Parameter Ambulatory Telemonitoring System Model CPC1S and how to connect it to a patient. Please follow the instructions, as they will assist you in getting the best possible results with the product. If you have any further questions, remarks for the supplier, or if you need more assistance, please contact the manufacturer or your local representative.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Contents

Intended Use and Usage Area	6
User Group	7
Contraindications:.....	7
Side Effects:	7
Complications:	7
Product life	8
Device Operating Instructions	8
1. Patient Preparation.....	8
2. Device Installation.....	8
ECG Electrode Patch Application Procedure	9
ECG Patch Configurations Available:.....	9
Important Notes:.....	9
Sensor Attachment.....	10
Logging into the “Check Point Cardio” telemonitoring system	10
Exiting the “Check Point Cardio” telemonitoring system.....	11
Choosing a medical center	11
Adding a new patient to the telemonitoring system	12
Adding a new visit to an existing patient in the telemonitoring system.....	14
Patient management.....	15
Monitoring session options.....	16
Patient data	16
Observation	17
Raw Data.....	18
Report review	19
Conclusions	20
Write conclusion	20

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Edit Patient File.....	21
Alarm settings.....	21
Medical Centers.....	25
Medical Staff.....	27
Accounts.....	29
General Alert Board	31
History of patient vitals.....	32
Disconnect.....	33
Questionnaires	33
Questions	35
Roles.....	37
MONITORING SCREENS.....	38
ICU monitoring.....	39
ICU monitoring screen for an individual patient	40
Observed parameters.....	41
Alerts	42
Other Possibilities.....	43
Observation dashboard.....	43
Individual patient observation	44
Alerts – Observation Dashboard	46
Additional “Observation Dashboard” options.....	46
Tracking	47
SYSTEM.....	49
Devices list	49
AUDIT LOGS.....	51
View log	51
Activating the Communication Device	53
Patient Monitoring	54
Automatic Alert and Patient Physiological Parameter Reporting System.....	55

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Alerts for Clinical Response According to the NEWS2 Standard:	56
Device Specifications	58
Warnings and Precautions.....	59
Liability Notice.....	60
Terms of Warranty.....	60
Device Insulation Diagram.....	60
Manufacturer's EMC Declaration	61
Symbols and Meanings in the Label Prepared in Accordance with EN 15223-1 Standard	63

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	



Product Name: Multiple Physiological Parameter Ambulatory
Telemonitoring System Model CPC1S

Company Name	Check Point R & D Ltd.
Authorized Person / Appellation	Ventseslav Kolev
Company Address	37 Osvobojdenie Str, 6100 Kazanlak / Bulgaria
Production Address	37 Osvobojdenie Str, 6100 Kazanlak / Bulgaria
Contact Number	+359888920108
Official Website / E-mail	https://www.checkpointcardio.com / info@checkpointcardio.com

Medical device model CPC1S



The system contains the following parts:

Medical device

- Main Sensor Unit -The primary medical device

Accessories:

- Battery Pack and backup battery - External battery used to charge the main unit
- 3 channel EKG patch- Patch for recording a 3-lead ECG signal

Non medical part:

- Adapter -Mains power adapter

The mobile application transmits data from the medical device BLE (Bluetooth Low Energy) to the file server via 3G/4G connection.

The power bank charges the internal battery of the device when connected to the second USB-C port.

The power bank must be connected only to a power source with 110–220 V, 50–60 Hz.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Intended Use and Usage Area

CPC1S is used to monitor vital physiological parameters (Electrocardiography, Non-invasive Blood Pressure, Respiratory rate, Pulse Rate and Body position and activity) of the human body by remote observation in continuous mode and real time and it helps the physician to diagnose the patient by following the vital physiological parameters. Medical Device CPC1S cannot replace medical examinations and consultations!

The device is designed for the observation of vital physiological parameters in real time. It is used for preventative follow-up, hospital follow-up, and after a hospital stay.

Tracking indicators – Use cases:

The indication of a CPC1S refers to the specific conditions or circumstances under which it is used. These systems are typically indicated for:

Chronic Disease Management: For patients with chronic conditions like heart disease, diabetes, hypertension, and respiratory disorders, continuous monitoring can help manage these diseases more effectively.

Post-operative Care: Monitoring patients after surgery to ensure that vital signs remain within normal ranges and to detect early signs of complications.

Elderly Care: For monitoring the elderly, especially those with health issues that require regular observation, to ensure timely medical intervention when necessary.

Remote Patient Monitoring: For patients who live in remote areas or have difficulty accessing regular healthcare services, providing a means to monitor their health condition remotely.

Research and Data Collection: In clinical research, collect data on the physiological parameters of participants over time.

In each of these cases, the telemonitoring system is used to gather and transmit health data in real-time or at specified intervals to healthcare providers, allowing for timely medical decisions and interventions.

User Group

The product is designed for patients of all sexes weighing over 10 kg, suitable for both sick patients and healthy individuals wishing to monitor their physiological status. The device is intended for diagnosis, early detection, and prediction of diseases and life-threatening conditions across all medical specialties. With the data provided, physicians can promptly adjust medication regimens and rehabilitation plans for various patient groups—including children and elderly patients— in both hospital and outpatient settings.

Contraindications:

- ⌞ The device is not classified as waterproof. Do not use the device in wet environment such as: rain, bathrooms, showers and other areas where water could reach the device. The device works normally at the perspired patients.
- ⌞ The device is not defibrillator protected. Do not use it with external (non-implanted) defibrillators. For safety, always remove electrodes, patient lead wires, and the device from the patient before defibrillation.
- ⌞ The device is not intended for use in infants weighing less than 10 kg.
- ⌞ The device is not capable of recording ECG signal in the presence of implanted pacemaker pulses, the function of the ME equipment can be adversely affected by the operation of an implanted pacemaker. So, do not use the device in presence of implanted pacemaker.
- ⌞ The device should not be used if the user has a burn or a wound on or near their wrist.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Side Effects:

- ↳ No known side effects have been reported.

Complications:

- ↳ After longer stay of the device it is possible the allergic reaction of the skin from ECG electrodes which are not part of model CPC1S medical device. The device must be reattached with new electrodes every 2 days to a different place on the skin

Product life

The life expectancy of the Multiple Physiological Parameter Ambulatory Telemonitoring System Model CPC1S unit is five (5) years.

1. Device Installation

- ↳ Ensure the device is powered up and running by attaching the battery pack to the main processing unit.
After attaching the battery pack make sure that the green light on the main processing unit (which is also the location of the SOS button) is blinking.

Attach the 3M Red Dot Electrodes to each Node on the CPC1S

NB: If the patient is male and has a hairy chest, depending on the amount of hair it can potentially lead to a lower quality ECG signal. To avoid this, shave the application areas on the patient's body with a razor to ensure optimal signal quality

Attach the medical device (see Pic. 1) to the patient's body.




(Pic.1)

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Logging into the “Check Point Cardio” telemonitoring system

Enter the access address for the “Check Point Cardio” telemonitoring system in your “Google Chrome” browser’s address line on your personal computer.


Choose the interface language for the “Check Point Cardio” telemonitoring system from the drop-down menu. Enter your “Username” and “Password” to access the system. Select “Enter”.



Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Exiting the “Check Point Cardio” telemonitoring system

To exit the system, select the “Exit” button.



Choosing a medical center


By choosing a medical center, the observer selects which patients from which medical center will be under observation. From the drop-down list choose the medical monitoring center.



Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Adding a new patient to the telemonitoring system


Select „Patient Admission Form “.



Fill in all the required data in the appropriate tabs (“Basic info”, “Initial risk assessment”, “Lab and microbiological data”, “Image Diagnostic”, “Diagnosis and medications”, “Confirm”). Navigate by pressing the “Next” or “Previous” button on each of the tabs or select the appropriate tab directly.

The patient session cannot be added to the system if:

- Any field marked with a “*” sign and red outline is not filled in.
- One or more of the fields marked with a “*” sign contains invalid data.






Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

When you've filled in all the mandatory fields, the "Save" button in the "Confirm" tab will become active.

Click the "Save" button to add the patient session to the system and agree to the GDPR declaration by clicking the "Accept" button.

If you do not accept the GDPR declaration, the patient cannot be added to the system.






Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	


Adding a new visit to an existing patient in the telemonitoring system

If additional patient observation is required (for patients who have already been added to the system once), a new visit should be created.

To do this, navigate to the previous patient session in the “Patient management” menu (the provided filters can be used for easier searching) and press the red “+” button (next to “Patient ID”)



To finish adding the new visit to the existing patient, the user needs to fill in all the required fields once more and click the “Save” button (in the “Confirm” tab).




The screenshot shows the 'Add new patient' form in the CHECK POINT Cardio system. The 'Confirm' tab is selected and highlighted with a red box. The form contains various fields for patient information, including 'Oxygen', 'Thermometer', 'Notes', 'Device MAC', 'Operation mode', 'Doctor's name', 'Phone number', and 'Email'. At the bottom, there are two orange status bars: '54 % High cardiac risk' and '63.15 % High general risk'. A blue bar at the bottom states 'Please fill all the required data in order to ensure proper medical care!'. A red box highlights the 'Save' button at the bottom right.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Patient management

This menu contains patients that are currently under observation, as well as patients that are no longer being monitored. To access it, select “Patient Management” from the sidebar menu.



In this menu, you can look for a patient manually or by applying one or more filters to narrow down your search. You can apply these filters by making your selection in the relevant box and pressing the “Apply Filters” button.

To see a list of all patients, active and inactive, press the “Show all” button.

To see a list of patients that are currently under observation, check the “Observed” checkbox and press the “Apply Filters” button.

For each patient, the following information is shown: “Patient ID” (generated automatically by the system), “Patient’s name”, “Birth date”, “Gender”, “Email”, “Phone”, “Address”, “City”, “Country”.

Status - This field shows whether a given patient is under observation



The patient is under observation. The color is determined by the system based on the calculations of the patient’s general and cardiac risk. Patients with a **low** general and cardiac risk will have their status denoted in **green**, patients with a **medium** general and/or cardiac risk will have their status denoted in **yellow** and patients with a **high** general and/or cardiac risk will have their status denoted in **red**.



The patient is currently not under observation.

To see additional options and information related to a specific patient, click on the blue arrow next to the patient.

	700700008102	Checkpoint cardio	Мария Трифонова Щерева	1970-07-13	Female	w@a	08 78 98 15 66	бул. 23-ти ПШП 65	Казанлык	Bulgaria	
ID	Start date	End date	Device Mac	Device Option	Operation Mode	Priority	Risk	Status			
11628	2025-03-27 17:04:14		20:06:23:93	3 lead ECG	Real time streaming	Low	Low				
Patient data Observation Raw data Report review Export Vit. Params Export position Conclusions Write conclusion Short report Edit Patient File Alerts/History History patient vitals Disconnect/Patient											

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	


Additional patient information includes:

- **ID** – Session number.
- **Start date** – The date and time when patient monitoring began.
- **End date** – The date and time when patient monitoring ended.
- **Device Mac** – The serial number of the device attached to the patient.
- **Device option** – The number of ECG leads (depends on the equipment used by the patient)
- **Operation Mode** – Tells you if the device is set to transmit in real-time or every two minutes (depends on the choice made by the user when registering the patient session)
- **Priority** – Automatically determined by the system, based on the patient's general and cardiac risk
- **Risk** – Automatically determined by the system, based on the overall health status of the patient at the time of admission
- **Status** – Shows whether the monitoring session is active or inactive.

Monitoring session options

Patient data

By choosing the "Patient data" option you will be taken to the following screen.



Using the tabs shown you can view the summarized basic information of the patient after they have been placed under observation and from here you can also assign a new: diet, medicinal therapy, activities and questions for the patient to receive, answer and perform.


- **Basic Info** tab – Here the user can view a summary of the patient's information as it was entered into the system during registration, as well as view all additional activities that have been assigned and whether the patient has completed these.
- **Notes-Operator** tab – Here the user (or operator) can add notes that can be seen by himself, as well as other op-

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

erators who have access to the system and this patient. Each note can be individually defined by the operator that creates it as “visible to others” or not.

- **Questionnaire-patient** – Here the user can check the questions assigned to the patient and the patient’s answer to these questions and if the patient has answered at all.
- **Measurements (CHARTS)** – Here the user can see a graphical representation of the added measurements.
- **Measurements (TABLES)** – Here the user can see all added measurements in a table format.
- **Clinical Plan** – Here the user can use the options provided in this tab to add or modify the patient’s: “Medication”, “Diet”, “Additional medical procedures”, “Activities” and “Questionnaires”
- **Calendar** – Here the user can view a broader image of the patient’s clinical plan in the form of a personalized calendar. Clicking on each calendar date will provide detailed information about the patient’s duties pertaining to the clinical plan for that specific timeframe.

To add measurements to the patient file the user can press the “Add measurement” button, select the desired measurement and fill in the required information. After this, the user must click on the “Save” button to finalize the process.



Add measurement

☐ Blood Pressure
 ☐ Blood sugar (mmol/L)
 ☐ Temperature (°C)

☐ Weight (kg)
 ☐ Height (cm)
 ☐ Respiration Rate


☐ Heart Rate (BPM)
 ☐ Oxygen saturation

Type measure	Value	Data
Save Close		

To export any of this data, simply select the “Export” button or select the print icon.

Observation

The “Observation” option leads to the “ICU monitoring” screen of the patient.



The screenshot shows the "Patient management" interface. On the left is a sidebar with navigation options: admin, PATIENTS, USERS, and MONITORING SCREENS. The main area displays a list of patients with columns for Patient ID, Hospital, Patient's name, Birth date, Gender, Email, Phone, Address, City, Country, and Status. Below this is a detailed view of a patient's data, including a table with columns for ID, Start date, End date, Device Mac, Device Option, Operation Mode, Priority, Risk, and Status. The "Observation" tab is highlighted in the patient's data view.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Raw Data

By selecting the “Raw data” option, the user will be taken to a new menu where they can track the history of the patient’s vital parameters and ECG strips.

Alternatively, you can navigate to the raw data menu by pressing the  icon in the “ICU monitoring” screen of an individual patient.




To navigate through the recorded data, you can set the date, time and the temporal range of the observation interval (15 or 60 minutes). To go forwards or backwards through the recorded data by your selected interval, use the corresponding arrows next to the “Range” option. To navigate using your set parameters, press the “Load” button.



Use the   buttons to label noise or pathologies.

From the list that appears after pressing the corresponding button, select the pathology or noise.



Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Use your mouse to scroll through the ECG graph and specify the start and end point of the pathology. Each ECG strip can be labelled with multiple pathologies or noises. If you'd like to delete an already labeled pathology or noise, select it with your mouse and press the "Delete" or "Backspace" key on your keyboard.

To add points, intervals or show vital parameter values on the ECG strips, use one or more of the following buttons:



To view additional ECG strip information, select the "Show All" button.



To save the ECG strip press the  button. This saves the ECG strip to the patient's conclusion. When this operation has been completed successfully, the button color will change: .

To export the data for the corresponding ECG strip (in .xlsx format) or print the selected ECG strip, use one or both of the following buttons:  .

To make more detailed analyses or leave a comment on an ECG strip, use the corresponding tools from the following menu:



Report review

The "Report Review" option allows the user to generate a graphical representation of the monitored vital parameters (in a .pdf format).

From the "Date" drop-down menu, select the date of interest for which you would like to receive a graphical representation of the patient's vital parameters. Afterwards, from the "Types" drop down menu, select one, several or all the vital parameters which of interest to you and finally press the "Review" button to generate the report.


Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

You can use the built-in .pdf instruments to zoom in and out of the report, rotate the graph(s), save the graph(s) in a .pdf format to your PC or print the graph(s).



Conclusions

The “Conclusions” option allows you to view all the patient’s created conclusions. Conclusions are presented in a .pdf format. To open a conclusion, simply select it from the menu.



Write conclusion

This option allows you to generate a conclusion for the period during which the patient was monitored.

- Diagnoses – This is where patient diagnoses are shown. The field is automatically filled with diagnoses that were specified when the patient was originally being added to the system for monitoring (patient admission) and can also be expanded on by the user if needed.
- Medication – This is where the current medication that the patient is taking is shown. Please note that this field is automatically filled with the medication that was specified when the patient was originally added to the system for monitoring (patient admission) and can also be expanded on by the user if needed.
- Use the provided checkboxes to specify which parameters you’d like to be included in the conclusion. Some of the parameters need to be manually selected from a drop-down menu.

SVE: ☒ under 5 per hour

VE: ☒ ▼

SVT: ☐ Lone class I

VT: ☐ Lone class II

AF/Af: ☐ Lone class III

AV blo: ☐ Lone class IVa

SA blo: ☐ Lone class IVb

ST: ☐ Lone class V

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

- Pathology – Select one or more pathologies from the drop-down menu, to be added to the conclusion. These selections will then be added to the “Conclusion” field.
- Conclusion - Write down the patient’s conclusion in this field.
- Recommendations – Any recommendations for the patient should be written here.

All recorded pathologies and ECG strips from the “Raw data” menu that were saved are then graphically represented in the conclusion. If you wish to remove one or more of them, select the trashcan icon in the upper right corner of the ECG strip.



A graphical presentation of one or more of the vital parameters can be added to the conclusion by selecting the checkbox of the vital parameter which the user wishes to add to the conclusion.

Automatic report

Select all ☐

Table News: ☐ News Total Score: ☐ Pulse Blood Pressure: ☒ Respiration: ☐ Oxygen Saturation: ☐ Temperature: ☐ Heart Rate: ☒ Extrasystoles: ☐ Pathologies: ☐ Tachyarrhythmias: ☐ Bradyarrhythmias: ☐ Heart Rate Variability Analysis: ☐

Stress Index: ☐ Shock Index: ☐ Centralization Index: ☐ General Deterioration Score (Risk): ☐ Blood Pressure: ☐ Mean Arterial Blood Pressure: ☐ Cardiac Parameters: ☐

To finish creating your conclusion, press the “Create Conclusion” button.


Edit Patient File

This option allows you to change or update any information related to the patient after they have already been added to the system for monitoring.

Alarm settings


This option allows you to change and set custom alert settings for each of the measured vital parameters (Respiration rate, Heart rate, , Blood Pressure, Temperature) for your medical center.

To access the alarm settings screen, you must first select the blue bell, which is visible next to the name of your medical center:




Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Once you have selected the blue bell you will see the following screen, which will allow you to begin making changes to the general values of the measured vital parameters for the whole medical center. These settings will apply to all patients registered to your medical center:



	Red (low)	Orange (low)	Yellow (low)	Yellow (high)	Orange (high)	Red (high)
Systolic	81.00	95.00	100.00	136.00	140.00	220.00
Respiration	10.00		16.00	20.00		21.00
Pulse	40.00		50.00	91.00	111.00	131.00
Temperature	35.00		36.00	38.10	39.10	
SpO2 Scale 1	91.00	93.00	94.00			
SpO2 Scale 2	83.00	85.00	87.00	93.00	95.00	97.00


Clicking on either the up or down arrows will allow you to raise and lower the value, at which the specific alarm will appear:



	Red (low)	Orange (low)	Yellow (low)	Yellow (high)	Orange (high)	Red (high)
Systolic	81.00	95.00	100.00	136.00	140.00	220.00
Respiration	10.00		16.00	20.00		21.00
Pulse	40.00		50.00	91.00	111.00	131.00
Temperature	35.00		36.00	38.10	39.10	
SpO2 Scale 1	91.00	93.00	94.00			
SpO2 Scale 2	83.00	85.00	87.00	93.00	95.00	97.00

Above each of the measured parameters you will notice un-editable text that says either “Red, Orange or Yellow” followed by either “low” or “high”:

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	




Alarm settings [X]

	Red (low)	Orange (low)	Yellow (low)	Yellow (high)	Orange (high)	Red (high)
Systolic	81.00	95.00	100.0	136.0	140.0	220.0
Respiration	10.00		16.00	20.00		21.00
Pulse	40.00		50.00	91.00	111.0	131.0
Temperature	35.00		36.00	38.10	39.10	
SpO2 Scale 1	91.00	93.00	94.00			
SpO2 Scale 2	83.00	85.00	87.00	93.00	95.00	97.00

CANCEL **OK**

Red, Orange or Yellow signify the urgency of the alert, with the red color denoting the highest urgency, whilst yellow denotes the lowest urgency, with orange being the middle ground. Low and high respectively address whether the parameter is lower than what is considered normal (e.g red low 91) or higher (e.g red high pulse 131).

If you're happy with the changes you've made, select "OK", otherwise select either "Cancel" or the "X" in the top-right corner:




Alarm settings [X]

	Red (low)	Orange (low)	Yellow (low)	Yellow (high)	Orange (high)	Red (high)
Systolic	81.00	95.00	100.0	136.0	140.0	220.0
Respiration	10.00		16.00	20.00		21.00
Pulse	40.00		50.00	91.00	111.0	131.0
Temperature	35.00		36.00	38.10	39.10	
SpO2 Scale 1	91.00	93.00	94.00			
SpO2 Scale 2	83.00	85.00	87.00	93.00	95.00	97.00


CANCEL **OK**

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

If you'd like to edit the personal settings for a specific patient to accommodate for differences in their vital parameters that are either considered normal or need to be watched more closely, you can access that patient's personal alarm settings by first going to "Patient Management":



From here select "Observation"



Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	


Once you've entered the "Observation" screen, select the blue button with the cogwheel on it to access the patient's personal alarm settings:



Selecting the blue cogwheel button will open the same alarm settings menu as the general medical center alarm menu. Editing the vital parameters there is done using the same process as in the medical center alarm settings.


Medical Centers

The system has the option to add, remove and edit any additional medical centers. To access the "Medical Center" menu, follow these steps:



Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Select the “Medical Centers” option from the main dashboard. Once you’ve done this, you will be taken to the following screen:




The screenshot shows the 'Medical centers' screen. On the left is a sidebar with navigation options: admin, PATIENTS, and USERS. The main area has a 'Filters' section at the top with input fields for Medical Center ID, Name, Short Name, Email, and Phone, followed by 'Apply filters' and 'Show all' buttons. Below is a table titled 'Medical centers' with a blue '+' icon in the top left corner. The table has columns: NR, Medical Center Name, Medical Center Short Name, Medical Center Email, Medical Center Phone, and Actions. The first row shows NR 2027, Name Bisturmed, Short Name Bisturmed, and Actions with 'Edit' and 'Delete' buttons highlighted by a red box. The second row shows NR 2026, Name DIGD - D2Oost, Short Name DIGD - D2Oost, and Actions with 'Edit' and 'Delete' buttons. The third row shows NR 2025, Name Oxford, Short Name Oxford, and Actions with 'Edit' and 'Delete' buttons.

From this screen, you can take the following actions:


- You can add a new medical center by pressing on the blue “+” icon
- You can edit an already existing medical center by selecting “Edit” under the “Actions” tab
- You can remove an already existing medical center by pressing on “Delete” under the “Actions” tab

If you choose to add a new medical center, you will be taken to the following screen



The screenshot shows the 'New Medical center' form. On the left is the same sidebar as the previous screen. The main area has a form with four input fields: Medical Center Name, Medical Center Short Name, Medical Center Phone, and Medical Center Email. The 'Medical Center Name' and 'Medical Center Short Name' fields are marked with a red asterisk (*). Below the form is a blue 'Save Medical Center' button.

In this menu you need to fill in at least the two mandatory fields marked by a red *. Afterwards, when you’re happy with the information you’ve put in, press “Save Medical Center” to create the new medical center.




This screenshot is identical to the previous one, showing the 'New Medical center' form. The 'Medical Center Name' and 'Medical Center Short Name' fields are highlighted with red boxes, and the 'Save Medical Center' button is also highlighted with a red box.


Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Medical Staff

The system has the option to add a list of medical staff to your medical center. To do so, select the “Medical Staff” option from the main dashboard:



After you’ve selected “Medical Staff” you will see the following screen:




From this screen you can:

- See already existing medical staff for your medical center.
- Add new medical staff using the blue + button.
- Edit an already existing medical staff member using the “Edit” button under the “Actions” tab.
- Delete an already existing medical staff member using the “Delete” button under the “Actions” tab.


Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

If you choose to add a new medical staff member by selecting the blue + button, you will see the following screen:



The screenshot shows the 'New Medical Staff' form. The left sidebar contains navigation links: admin (Project Administrator, Profile, Exit), PATIENTS (Patient admission form, Patient management, Questionnaires, Questions), USERS (Medical centers, Medical staff, Accounts, Roles), and MONITORING SCREENS. The main form area has the following fields: Email *, First name *, Surname, Family name *, Phone, Medical staff type *, Academic Title, and Experience *. The 'Save Medical Staff' button is at the bottom left. Red boxes highlight the mandatory fields: Email, First name, Family name, Medical staff type, Experience, and the Save Medical Staff button.

From here you need to fill in all the mandatory fields marked by a red *, as well as select an option from the mandatory drop-down menus:




This screenshot shows the 'New Medical Staff' form with the 'Medical staff type' dropdown menu open. The dropdown menu lists two options: 'Doctor' and 'Nurse'. The mandatory fields (Email, First name, Family name, Experience) and the 'Save Medical Staff' button are still visible. The 'Medical staff type' field is highlighted with a blue border.

Once you have filled all the mandatory fields, select the blue “Save Medical Staff” button to save your new staff member.


Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Accounts

The new system allows the creation of multiple accounts able to access the same medical center/s. To view, edit and create new accounts for your medical center, select the “Accounts” option from the main dashboard:



Once you’ve selected the “Accounts” tab you will see the following screen:




NR	User Name	User Type	Medical Center	First name	Family name	Status	Actions
607	PAAppleTest	admin	Project Administrator	Apple	Test	☒	
602	KatiePA	admin	Project Administrator	Katie	Weller	☒	
512	AstenCante	admin	Project Administrator	Asten	Cante	☐	

From this screen you can:

- View already existing users, their username, the type of user, the medical center they are registered to, their status, which can be either active or inactive (marked with a blue checkmark for active and an empty box for inactive).
- Edit an already existing user by pressing on the blue paper and pen button to the right of the user you’d like to edit under the “Actions” tab.
- Delete an already existing user by pressing on the blue trashcan button to the right of the you’d like to delete under the “Actions” tab.
- Add a new user by pressing on the blue “New User” button.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

If you choose to add a new user, you will be taken to the following registration form:



In this registration form you need to perform the following steps to successfully create a new user:

- ┆ Set the status of the user to “Active”
- ┆ Select the user type from the drop-down menu (It is recommended you create an admin account before you create a medical staff account).
 - ┆ Choose the roles you would like this user to have. These roles are best described as the rights to certain menus and features that the user you are creating will have (e.g. adding/removing patients, accessing the “Accounts” tab, etc.)
- ┆ Select the medical centers you would like this user to be able to view.
- ┆ Choose and enter a username for the new user.
- ┆ Choose and enter a password for the new user, which complies with the requirements listed below the password field. If you’d like the system to generate a confirmed and secure password for you, click on the grey “Generate” button next to the passwords field.
- ┆ Confirm the password that was generated or chosen.
- ┆ Enter the user’s first name.
- ┆ Enter the user’s surname (optional).
- ┆ Enter the user’s family name.
- ┆ Enter the user’s email address.
- ┆ Enter the user’s phone number (optional).

Once you’ve completed all the above steps correctly, click on the blue “Save” button and your user will be created.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

General Alert Board

This option allows the user to observe all alerts generated by the system, in real-time, for their selected patient or all active patients.



The screenshot shows the 'GENERAL ALERT BOARD TEST TEST, 04' interface. It is divided into two main sections: 'NEWS ALERTS' and 'CARDIAC ALERTS'.

NEWS ALERTS:

ID	ALERT TIME	EVENT	FOLLOW-UP ACTION	RESP. TIME
70689	23.02.12:17:54	Respiration rate value is high	[Dropdown]	
70688	23.02.12:12:11	Heart rate value is high	[Dropdown]	
70687	23.02.12:10:18	Low device battery	[Dropdown]	


CARDIAC ALERTS:

ID	ALERT TIME	EVENT	FOLLOW-UP ACTION	RESP. TIME
70709	23.02.12:07:52	Atrial couplets	[Dropdown]	
70708	23.02.12:07:56	Atrial couplets	[Dropdown]	
70707	23.02.12:07:58	Atrial couplets	[Dropdown]	
70706	23.02.12:08:15	Atrial couplets	[Dropdown]	
70705	23.02.12:08:15	Atrial extrasystoles in trigeminy	[Dropdown]	

Both “NEWS alerts” and “Cardiac alerts” are displayed for the current user.

Alerts in the “NEWS alerts” section are raised when some of the monitored parameters exceed preset limits. To specify the limit for when the alert will be generated for one or more parameters press the  button.

You can also divide alerts by color based on the preset NEWS 2 limit.



The 'ALARM SETTINGS' dialog box allows users to configure alert thresholds for various parameters. The settings are organized by color-coded severity levels: Red (low), Orange (low), Yellow (low), Yellow (high), Orange (high), and Red (high).

Parameter	Red (low)	Orange (low)	Yellow (low)	Yellow (high)	Orange (high)	Red (high)
Systolic BP	90.00	100.00	110.00			220.00
Respiration	8.00		11.00		21.00	25.00
Pulse	40.00		50.00	91.00	111.00	131.00
Temperature	35.00		36.00	38.10	39.10	
SpO2	91.00	93.00	95.00			

Buttons: CANCEL, OK

Please note that alerts related to problems with patient measuring equipment are also raised in this section.

Alerts in the “Cardiac alerts” section are sent when a pathology is labeled in the “Raw data” screen of the current patient. Based on the type of pathology, the alert is marked with a different color.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

When the alert is sent (in both sections) it is expected that the user will evaluate the alert and take necessary action. These are predefined in the “Follow-up action” drop-down menu.

ID	ALERT TIME	EVENT	FOLLOW-UP ACTION	RESP. TIME
70689	23.02.12:17:54	Respiration rate value is high	Continue routine observation	24.02.11:34:16
70688	23.02.12:12:11	Heart rate value is High	Urgent assessment from the ward	24.02.11:34:56
70687	23.02.12:10:18	Low device battery	Continue routine observation Urgent assessment from the ward doctor Emergent assessment from the clinical team Reject	

The date and time when the follow-up action is taken is recorded in the “Resp. time” section. Please note that when a “Follow-up action” is taken for specific alert, that alert is deleted from the “General alert board”.


By using the  button, the user can see all rejected alerts.

History of patient vitals

This observation screen allows the user to track the values of the patient’s vital parameters in real-time, as well as simultaneously view the history of their vital parameters for a specific date and time.

The measured or displayed parameters depend on the type of CPC device the patient is currently wearing.

To see the recorded values of the measured parameters, select the date from the list. Use the “Previous” or “Next” buttons to change the date. For a more detailed digital representation of the measured parameters, highlight one of the graphs with your mouse cursor and use the mouse to further navigate in the graph.




Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	


Disconnect

This option allows the user to disconnect the patient's session from the system, when the patient is no longer during observation.

To do this, select the "Disconnect" option. Afterwards, select the "Accept" button to confirm you wish to disconnect the patient from the monitoring system or "Reject" if you wish to cancel the disconnection process.



Finally, close the information window by pressing "OK".




Disconnected patients can only be re-added to the system by creating a new patient session.

Questionnaires


This option allows the user to create custom questionnaires or use preset ones, which can then be assigned to individual patients or a group of patients.

From the main dashboard, select "Questionnaires".




Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

The user can use the options provided to filter through the list of already created questionnaires. Use the “Show All” button to display all the available questionnaires.



To create a new questionnaire, select the **New questionnaire** button.



- Questionnaire name – Assign a name to your questionnaire.
- Type – Specifies the questionnaire type: “Group” or “Personal”.

Group questionnaires will be assigned to all patients who belong to a predefined disease specified in the “Patient group” and “Patient subgroup” drop down lists in the “Diagnosis and medications” tab when the patient was originally being added to the system.

Personal questionnaires will only be assigned to individual patients whose name or anonymized callsign can be selected via a drop-down menu.

- Notification settings** – This allows for notification intervals and settings to be adjusted further by the user. The following options are available: “No notification”, “Weekly”, “Days of the month” (chosen by the user), “In a few days” and “A certain date”. Please note that if notifications are enabled, the patient will be notified of each questionnaire on their mobile device.
- Add questions** – Here the user can add questions to the questionnaire by selecting them from a drop-down menu or they can create new questions using the “New Question” button.

Questions created through the “New Question” button are saved and can be viewed in the “Questions” menu.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	


To save a newly created questionnaire, select the “Save” button. To

edit an already existing questionnaire, use the  button.


To delete an already existing questionnaire, use the  button.

Questions

This option is used to create individual questions which can then be added to a questionnaire. From the main dashboard, select “Questions”.




The user can use the options provided to filter through the list of already created questions. Use the “Show All” button to display all the available questions.



The screenshot shows the 'Questions' page. It features a 'Filters' section with a text input field labeled 'Question'. Below the input field are two buttons: 'Apply filters' and 'Show all'.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

To create a new question, select the **New question** button.



- **Write the question** – In this field you must write down the question. When the question is assigned to the patient or patient group through “questionnaires”, it will be sent to the patient. The action that the patient must take to answer the question is dependent on the question’s content or type.
- **Content of the question** – From here you must select the content of the question:
- **Free text field** – The patient will have to answer the question using free text, before sending their answer back using the “Send” button.
- **Several answers to choose** – This option allows for several multiple-choice questions and answers to be created. For each created choice the “Type of answer” needs to be specified. The following types of answers are available: “Positive”, “Negative” and “Undefined”.

The patient receiving this type of question will have to choose one of the available answers to send back.

- **Measure parameter** – When the patient receives this type of question, they are expected to measure the specified vital parameter using the appropriate measuring equipment and then fill in the value of their measurement before sending the answer back to the system.

Patient answers can be seen in the “Questionnaire-patient” tab. To

save newly created questions, select the “Save” button.

To edit already existing questions, press the following button  To

delete already existing questions, press the following button 



Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Roles

This option is used to manage account roles, which allow or do not allow the user to perform specific actions within their assigned medical center.


Each account or “user” in the system has one or more roles attached to them. In turn, each role has permissions assigned to it. These permissions determine the amount of access the user is allowed to have within the system.

Select “Roles” from the dashboard to access the “Roles” menu.



Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

To create new role, select the “New role” button.



- Name – This field requires that you assign a name to your new role.
- Description – This field allows you to attach a short description the role you are creating.
- Assigned Permissions – All available role permissions are listed here. To add these permissions to your role, press the  button. To remove any permissions, press the  button.

To finish creating the role, select the “Save” button.

To edit an already existing role, press the  button. To delete an already existing role press the  button.

MONITORING SCREENS


The “Check Point Cardio” monitoring system offers different screen layouts for patient monitoring. Each one of the monitoring screens offers a different form of presentation of the monitored parameters. The user can choose whether to monitor multiple patients or see a more detailed view of an individual patient.

Information displayed on the monitoring screens depends on the used patient measuring equipment (device).

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

New ICU monitoring screen

To access the new ICU screen, press on the “New Screen” button in the leftmost toolbar:



Pressing on “New Screen” allows you to access the updated ICU view, from which you can observe active patients, expand their observation for more information, add new patients and switch through medical patients:



Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

To expand the patient view simply press on the patient you'd like to observe from the list to get a more detailed view of their trends, vital parameters, live ECG waveforms and more:




From this expanded patient view you can see the patient's live vital parameters, body position, NEWS score and their trend, raw, live and alert data and also technical information like the battery level of the communication device and the medical device. To switch through these various menus simply click on them. Furthermore, you can also scroll through the days during which the patient was under observation by selecting the appropriate date (dd-mm format) and review their past data.




Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

By accessing the “Realtime Waveforms” tab you can see the live ECG data of the patient. Depending on the type of device the patient is wearing, the number of ECG boxes changes.




From the “Raw Data” tab, you can access the raw data of the patient. In this tab, you can scroll through the raw data using the “next” and “previous” buttons or through the “Date” drop-down menu. Additionally, using the “Measurement” button, you can make take any necessary measurements of the patient’s ECG data directly in this tab. From here you can also download important ECG strips or print them out using the “Print” and “Download”




Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Through the “Alerts” tab you can view all automatic and manual alerts that are generated for this particular patient session with their appropriate timestamps. From here you can also set the personal alert threshold of a patient using the “Alerts Threshold Settings” button.



Date	Time	RESP (brpm)	SpO2 (%)	HR (bpm)	TEMP (°C)	cNiBP (mmHg)	Alert
19.08	10:36:00						Systolic BP is Low
19.08	10:32:00						Heart rate value is High
19.08	10:18:20						Respiration rate value is low
19.08	09:56:40						Respiration rate value is low
19.08	09:21:40						Heart rate value is High
19.08	09:06:00						Respiration rate value is low
19.08	08:42:00						Systolic BP is Low

Clicking on the “Add New Patient” button will open the short patient admission form where you can enter the data of a new patient and register them through the New ICU screen directly. To add the patient simply fill-in all the mandatory fields with the patient’s information and press “Save”.



Short patient admission

External patient ID:

Oxygen *: ☐ Yes ☐ No Thermometer *: ☐ Yes ☐ No

First name *:

Patient location:

Surname:

Device serial number *:

Family name *:

Urgent phone number *:

Birth date *: Gender *:

Height (cm) *: Weight (kg) *:


Systolic BP *: Diastolic BP *:

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

If your account has been given access to more than one medical center, you can scroll through the available options by clicking the white arrow next to the current medical center you are on and selecting the medical center you'd like to switch to from the drop-down menu.

12 ДКЦ	AmericaForBulgaria Devin	AmericaForBulgaria
AmericaForBulgaria Lovech	AmericaForBulgaria Vratsa	Kazanlak
ASPIRONIX	Asten Sante A Domicile	Asker Demo
Bisturmed	Bulgaria-TEX	Avant Maris Medical
Cigalah Gulf Medical Co	Comac Medical	Checkpoint cardio
Dedine 2	Dedine 2	Control - Communication
Demo2022		Center - Vraca
Fakultni nemcnice Olomuc	General University Hospital - Prague	DEMO
Home Care Center - Biala Slatina	Home Care Center - Krivodol	Home Care Center -
Home Care Center - Oriahovo	Home Care Center - Vidin	Belogradchik
HOME-ART	Human Byte	Home Care Center -
Institute for Clinical and Experimental Medicine Kepler Universitätsklinikum	IRHIS	Montana
LIDL	KraliMarko	Home Care Center - Vraca
	Longevia Health	IKVB
	Medica	
	Medtech	Loramed
Moja Klinika	Oborishte	Medintech
OneMed-DEMO	OstraMostra	MOBAL D-r Stefan
Petar Medical Center	Project Administrator	Cherkezov AD
Robin	Royal Liverpool Hospital	ONCOLOGY UKCLJ
SatHealth	Seha Virtual Hospital	Oxford
SimpleComply	Telecom Slovenia	Q16
		Sanjiivv
		Sheffield
		TSTest


By clicking on the “Calibrate” button when you are in the expanded patient menu you can adjust the initial BP value entered during patient registration by simply entering the new value in the Systolic and Diastolic BP fields and pressing “Save”.




Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Tracking

This menu is used to monitor the physical position of all currently active patients on “Google Maps”.



To access this menu, select “Tracking” in the main dashboard.




Based on the patients’ health risk, each positional marker is presented in the according color.

Use the provided menu (in the upper right corner) to filter the patients based on their health risk. Only the options that are checked will be displayed on the map.


Fit Bounds ☒ All ☒ ☒ Mobile team ☒ ☒ Very high Risk ☒ ☒ High Risk ☒ ☒ Middle Risk ☒ ☒ Low Risk

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

To receive additional information about the patient, click on their position marker.



Use this icon (in the upper left corner) to quickly switch between the following screens: “ICU monitoring”, “Observation dashboard”, “General alert board” and “Tracking”




Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

SYSTEM


Devices list

This menu is used to check which medical devices are assigned to a medical center. To access this menu, select "Device List" from the main dashboard.



The screenshot shows the main dashboard with a sidebar on the left. The sidebar has a 'MONITORING SCREENS' section with options like 'ICU monitoring', 'High risk patients', etc., and a 'SYSTEM' section with options like 'Sim Cards List', 'Devices list' (highlighted with a red box), 'Device Payment Info', 'Translations', and 'Phone app versions'. The main area shows a 'Dashboard' header with a 'General Statistics' chart and three summary cards: '11437 Total patients', 'N/A Conclusions per month', and 'Calendar Scheduling'.

When first accessing the menu, you will see all the currently available devices for your medical center.




The screenshot shows the 'Device list' menu. It has a 'Filters' section with input fields for 'Device Mac' and 'Device Patch Number', and a 'Show all' button. Below is a table of devices with columns: NR, Device Mac, Device Patch Number, SIM Tel, SIM ID, Device model, Medical Center, Device Status, Medical staff, SPO sensor type, Registration date, Manufactured date, Activity, and Actions. The table contains 6 rows of device data.

NR	Device Mac	Device Patch Number	SIM Tel	SIM ID	Device model	Medical Center	Device Status	Medical staff	SPO sensor type	Registration date	Manufactured date	Activity	Actions
807	201113494	34			CPC3	Demo2022	In warehouse	Dimitar RADEW	Finger sensor	2025-03-21	2025-03-21		Edit, Delete
806	221113492	34			CPC3	Demo2022	In warehouse	Dimitar RADEW	Finger sensor	2025-03-21	2025-03-21		Edit, Delete
805	20092910	26			Nightingale CPL125	Karchenska Institut	to client	Karolina KAROLINSKA	Ear clip sensor	2025-03-20	2025-03-20		Edit, Delete
804	25.02.0172	01			CPC 25	Project: Administrator	to client	Evgenia Cerepuchova	Finger sensor	2025-03-12	2025-03-12		Edit, Delete
803	25.02.0171	01			CPC 25	Project: Administrator	to client	Evgenia Cerepuchova	Finger sensor	2025-03-12	2025-03-12		Edit, Delete
802	25.02.0170	01			CPC 25	Project: Administrator	to client	Evgenia Cerepuchova	Finger sensor	2025-03-12	2025-03-12		Edit, Delete


Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

To search for specific devices, use the provided filter options. To see all devices registered to the medical center, select the “Show all” button.



To add a new device to your medical center, select the blue “+” button.

After selecting the blue “+” button you need to follow these steps to finish adding the new device to your medical center



- └ **Device Mac** – Enter the serial number of the device. The serial number or “MAC” is the unique physical address of a device. Please follow the specified MAC address format when filling in the field.
- └ **Device Batch Number** – The field is related to the device batch number and is filled-in automatically by the system, based on the MAC address of the device.
- └ **Medical Center** – Select the name of the medical center that the device will be assigned to from the drop-down.
- └ **Device model** – Select the device model from the drop-down menu.
- └ **Medical staff**– Select the medical personnel that will have the device assigned to them from the drop-down menu.
- └ **sensor type** – Select the type of sensor that will be used with this device.
- └ **Device Status** – Select the status of the device:
 - **to client** – This status indicates that the device is currently in use by a patient. Devices that have this status attached to them cannot be selected from the “Device” drop-down menu in the patient admission form.
 - **in warehouse** – This status indicates that the device is currently not in use and can be selected from the “De- vice” drop-down menu in the patient admission form.
 - **in service** – This status indicates that the device is currently not available because it is undergoing maintenance at our service center. It will not be available for selection in the “Device” drop-down menu in the patient admis- sion form.

Please note that the “Device status” is automatically changed from “in warehouse” to “to client” when the device is being used by a patient. When the patient is disconnected, the device status will automatically change from “to client” to “in warehouse”.

- └ **Manufactured date** – Specify the manufacturing date of the device in this field.
- └ **Registration date** – Specify the date that the device was added to the system for the first time.

To edit already existing device information, select the Edit button. After you’re finished editing, select the “Save” but- ton.

To remove a device from the medical center, select the Delete button.


Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

AUDIT LOGS

View log


The “View Logs” menu is used to track user activity in the system for a specific medical center.

To access it, select the “View Logs” option from the main dashboard.



All the system information is stored in different database tables. Any action which the user takes regarding patient information such as adding, deleting or updating it anywhere in the system, is tracked.

You can use the provided filters to refine your search results.



The screenshot shows the 'Logs' page. At the top is a header with a menu icon and the word 'Logs'. Below the header is a 'Filters' section with four dropdown menus: 'Date', 'Action', 'Table', and 'Users'. To the right of these menus are two buttons: 'Apply Filters' and 'Show all'.

- └ **Date** – Set the date you are interested in from the drop-down menu.
- └ **Action** – Set the type of actions you are interested in from the drop-down menu. The following user interactions can be tracked:
 - **insert** – These types of actions are related to inserting information in database tables.
 - **update** – These types of actions are related to updating any existing information in specific database tables.
 - **delete** – These types of actions are related to deleting existing information in specific database tables.
- **Table** – Select the table for which you want to track changes in the drop-down menu.
- **Users** – Select a specific user or users whose actions you wish to track from the drop-down menu.

Once you’re done with your filter selection, select the “Apply Filters” option to apply them.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

To receive more detailed information regarding an action, select it. Information for the “Old value” and “New value” change is then displayed.

Logs						Export	Print
Date	User	IP	Description	Operation	Table		
2025-03-28 12:32		94.26.9.235	Update record of table "VitalAnalysisAlertFirebase"	update	VitalAnalysisAlertFirebase		
Old value		New value					
VitalAnalysisAlertsSent: 0		VitalAnalysisAlertsSent: 1					
2025-03-28 12:32		94.26.9.235	Update record of table "VitalAnalysisAlertFirebase"	update	VitalAnalysisAlertFirebase		
2025-03-28 12:31		94.26.9.235	Update record of table "VitalAnalysisAlertFirebase"	update	VitalAnalysisAlertFirebase		

The following information is available for each log:

- **Date** – Date and time when the parameter was changed.
- **User** – System user who made the change.
- **IP** – The IP address from which the user has accessed the system and made the change.
- **Description** – Specifies which table and what kind of action was executed.
- **Operation** – Specifies what kind of action was executed on the parameter by the user.
- **Table** – Specifies which tables were affected by the change.

To export or print the information in this log, select the “Export” or “Print” button respectively.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Activating the Communication Device

Switch the communication device on by holding down its power button. The communication device will then start, and you will be taken to the following screen.




This message confirms that the device has successfully connected to the server.



If yellow alert messages appear on the main screen, press the **Alert** button to view details about the issue.



If the device cannot connect to the server, swipe down on the screen and verify the WiFi or 3G connection status, as shown in the image above.



The communication device can switch between WiFi and 3G networks as needed.



The device temporarily stores data packets in a local database. The buffer value indicates how many data packets have not yet been sent to the server. Once the connection to the server is restored, the buffer value will reset to zero.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

THE MESSAGES SHOWN ON THE COMMUNICATION DEVICE

"Device is not assigned"	The device is not receiving acknowledgment signals confirming that the data has been successfully delivered to the file server.
"ECG sensor is detached"	The ECG signals are of unacceptable quality. Possible reasons include detached sensors or expired electrodes.
"Device is detached"	The device is not attached to the patient. It may be placed on a table or in a box. In this case, vital sign values will not be displayed on the screen.
"Please charge the communication device"	The battery level of the communication device is below 15%.
"Please charge the medical device"	The battery level of the medical device is below 15%
"Connecting to server"	The device is transmitting data packets to the server, but acknowledgment signals confirming receipt have not been received.
"Connected to server"	The device has received acknowledgment signals from the server, confirming that the transmitted data packets have been successfully accepted.
"Device is not connected"	There is no connection between the medical device and communication device. A problem with Bluetooth pairing is possible. Press the lock/unlock button 10 times to check
"Device connecting"	The communication device is attempting to establish a connection with the medical device.
"Low quality Bluetooth signal"	The distance between the medical device and the communication device is too large, preventing the normal transfer of both historical (buffered) data and real-time data streams.
"Your heart is online !"	If this message appears, the device is correctly installed, and communication has been successfully established.

Patient Monitoring

The medical device transmits data to the communication device via a Bluetooth Low Energy (BLE) connection. Real-time values of vital parameters are displayed on the 5.8-inch screen, along with statistical information over selected time periods. Raw data from all measured signals is transferred to the file server, where it is stored, displayed, and available for further analysis.

After placing the patch and connecting the device, open the **ICU Monitoring** page to verify that all ECG leads are correctly attached and functioning properly. Once proper connectivity is confirmed, the communication device will display the message: **"Your heart is online."** To monitor all vital parameters, navigate to the main menu and select the patient's name.



Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

All valid ECG Leads are used to detect the patient's **heart rate**. QRS complexes and their corresponding R-wave peaks (Rmax), along with the timestamps of maximum QRS amplitudes, are identified using a specialized velocity method. This approach involves generating a special velocity (SV) signal for QRS slope detection and employing a set of adaptive beat-detection thresholds.

The Hilbert Transform method has a unique property of producing zero-crossings at peak-value locations. Using these zero-crossing points, corresponding R-peak locations are identified to detect pauses (**Pause Detection**).

The medical device model **CPC1S** supports the detection of **ST elevation**. The device does not measure atrial fibrillation (**AF**) or ventricular fibrillation (**VF**)

Automatic Alert and Patient Physiological Parameter Reporting System

The system integrated into the CPC desktop application for monitoring physiological parameters is based exclusively on the standardized early warning decision support framework known as **NEWS2** (National Early Warning Score 2).

NEWS2 is the latest version of the National Early Warning Score (**NEWS**), initially released in 2012 and updated in December 2017. It provides a standardized approach for assessing and responding to acute clinical deterioration.

NEWS2 has received formal endorsement from NHS England and NHS Improvement, establishing it as the primary early-warning system for identifying acutely ill patients—including those with sepsis—in hospitals across England.

5.1. The NEWS Clinical Observation Chart

To ensure a standardized and consistent approach to recording vital signs data, a color-coded clinical chart (**NEWS Chart**) has been developed and widely implemented across the NHS. This chart is designed for recording routine clinical observations and monitoring a patient's clinical condition. The primary objective of this chart is to promptly alert clinical teams to any significant clinical deterioration and to facilitate effective monitoring of patient recovery. The NEWS score assists in determining the urgency and scale of the clinical response required.

5.2 NEWS 2 Observation screen layout:

Chart 1: The NEWS scoring system

Physiological parameter	3	2	1	Score 0	1	2	3
Respiration rate (per minute)	≤8		9–11	12–20		21–24	≥25
SpO ₂ Scale 1 (%)	≤91	92–93	94–95	≥96			
SpO ₂ Scale 2 (%)	≤83	84–85	86–87	88–92 ≥93 on air	93–94 on oxygen	95–96 on oxygen	≥97 on oxygen
Air or oxygen?		Oxygen		Air			
Systolic blood pressure (mmHg)	≤90	91–100	101–110	111–219			≥220
Pulse (per minute)	≤40		41–50	51–90	91–110	111–130	≥131
Consciousness				Alert			CVPU
Temperature (°C)	≤35.0		35.1–36.0	36.1–38.0	38.1–39.0	≥39.1	

5.3. Definition of alerts and thresholds:

A score is assigned to each measured parameter, with the magnitude of the score corresponding to the degree of deviation from the normal range. These individual scores are then aggregated into a total score. For patients requiring supplemental oxygen to maintain recommended oxygen saturation levels, the total score is increased by 2 points. Alerts and their corresponding values are displayed according to the NEWS standard in green (NEWS 0), yellow (NEWS 1), orange (NEWS 2), and red (NEWS 3 or higher).

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

See the figure below for details:

Chart 2: NEWS thresholds and triggers

NEWS score	Clinical risk	Response
Aggregate score 0–4	Low	Ward-based response
Red score Score of 3 in any individual parameter	Low–medium	Urgent ward-based response*
Aggregate score 5–6	Medium	Key threshold for urgent response*
Aggregate score 7 or more	High	Urgent or emergency response**

* Response by a clinician or team with competence in the assessment and treatment of acutely ill patients and in recognising when the escalation of care to a critical care team is appropriate.

**The response team must also include staff with critical care skills, including airway management.

Alerts for Clinical Response According to the NEWS2 Standard:

Chart 4: Clinical response to the NEWS trigger thresholds

NEWS score	Frequency of monitoring	Clinical response
0	Minimum 12 hourly	Continue routine NEWS monitoring
Total 1–4	Minimum 4–6 hourly	- Inform registered nurse, who must assess the patient - Registered nurse decides whether increased frequency of monitoring and/or escalation of care is required
3 in single parameter	Minimum 1 hourly	Registered nurse to inform medical team caring for the patient, who will review and decide whether escalation of care is necessary
Total 5 or more Urgent response threshold	Minimum 1 hourly	- Registered nurse to immediately inform the medical team caring for the patient - Registered nurse to request urgent assessment by a clinician or team with core competencies in the care of acutely ill patients - Provide clinical care in an environment with monitoring facilities
Total 7 or more Emergency response threshold	Continuous monitoring of vital signs	- Registered nurse to immediately inform the medical team caring for the patient – this should be at least at specialist registrar level - Emergency assessment by a team with critical care competencies, including practitioner(s) with advanced airway management skills - Consider transfer of care to a level 2 or 3 clinical care facility, ie higher-dependency unit or ICU - Clinical care in an environment with monitoring facilities

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

1. Physical Device Removal

To remove the CPC1S from the patient, simply detach the device from the ECG electrodes, one by one.

⚠ Important: Do not pull directly on the CPC1S, as this can significantly shorten its lifespan.

2. Charging the Device

The CPC1S medical device is powered by a rechargeable Li-Ion battery pack.

To charge the battery pack:

- Connect the battery pack to the charging adapter and attach it to a power socket.

The LED indicators will begin flashing slowly from one LED, progressively lighting up until all LEDs are steadily illuminated, indicating the battery is fully charged.

3. Device Maintenance and Cleaning

It is recommended to clean the device and patch before each patient examination:

- Use a soft cotton cloth dampened with 70% alcohol solution.
- Wipe gently without applying excessive pressure, giving particular attention to any embossed or raised areas.
- When cleaning the device's side panel, ensure liquids do not enter the cable or connector outlet.

The CPC1S may be cleaned using either hydrogen peroxide (3%) or isopropyl alcohol (70%).

4. Proper Disposal of the Product

The CPC1S should be disposed of according to the institution's medical waste disposal procedure. It does not require a special disposal process.

5. Troubleshooting

If the communication device displays the message **"Device disconnected"** and the red IR sensor on the medical device is on, verify proper pairing:

- Press the **Lock/Unlock** button **10 times** consecutively.
- A list of nearby Bluetooth devices will appear. The paired medical device will be indicated with a star (★) next to its name.

6. Product storage and Operational Conditions

- Environmental Conditions:
 - Operating Temperature: 5°C~40°C Storage Temperature: -20°C~55°C
 - Humidity: ≤80%, no condensation in operation
 - ≤93%, no condensation in storage

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Device Specifications

MEDICAL DEVICE MODEL CPC1S	
Dimensions	42mm x 49mm x 16mm
Weight:	0.030 kg
Operation time:	24 h
Power:	Battery Pack 3.7 V 650 mAh
Storage memory	1 day
Power:	Battery Pack 3.7 V
Power check button:	No
ECG	
Leads	3-Channel
ECG	4 Electrodes
Sample Rate	250 Hz
Lead Detection	Yes
Respiration rate	Yes
Blood Pressure	60 – 250
Systolic:	40 – 150
Diastolic:	

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Warnings and Precautions

	Model CPC1S is a device for monitoring the vital parameters of the human body. Its main purpose is to help medical specialists to monitor their patients. The CPC1S medical device cannot replace medical examinations and consultations!
	The device is not classified as waterproof. Do not use the device in wet environments such as: rain, bathrooms, showers, and other areas where water could reach the device. The device works normally with perspired patients.
	The Device is not defibrillator-protected, do not use with not-implanted defibrillator. For safety reasons — Remove the electrodes, patient lead wires, and main box from the patient before defibrillation.
	The device is not intended for infants weighing less than 10 kg
	Do not let your medical device model CPC1S get wet – liquids could cause serious damage.
	Do not store medical device model CPC1S in a dirty or dusty place to prevent damage.
	Avoid dropping, throwing, or crashing your device into other objects.
	The device could be damaged if it is exposed to magnetic fields.
	Do not expose your device to high or low temperatures. This could shorten the life of the battery or melt some plastic elements. When the environment in which the device is located has rapidly changing temperatures, it can get wet and could damage the electronic circuit boards.
	Do not use corrosive materials or cleaners to clean the device.
	Do not paint the device.
	Do not open the device on your own and do not let unqualified people do so. The device contains small parts that could be broken, as well as the device itself.
	Always carry medical device model CPC1S in its bag to ensure proper functioning. The bag is designed to provide the best protection for your device.
	There are no user-serviceable components inside the system. Authorized service personnel should do all internal troubleshooting and repair or replacement using only parts and accessories approved by Check Point R&D.
	Periodically check all connector cables for damage. Do not operate the equipment if the integrity of these components is questioned.
	Only people who are qualified, competent, and have the right knowledge about laws regarding telemedicine should operate this device.
	The device is not intended to be used with HF surgical equipment
	The Bluetooth® transmission may be influenced in different environments by materials such as walls, metallic doors, steel wire netting and an MRI or a CT environment.
	Never place the battery packs of the CPC1S in heating devices such microwaves and ovens or electric heaters. Batteries can explode if heated. Do not break the integrity of the battery pack. Be careful not to expose the battery to external pressure - this could cause short circuiting and overheating.
	The battery pack may heat while charging.
	Proximity to magnets, metal detectors, high-voltage electrical wires, and electrical appliances such as shavers, toothbrushes, and hair dryers may have an influence on the signal quality.
	Excessive perspiration may cause the electrodes to loosen or detach completely.
	Inaccurate measurements and readings may be caused by autoclaving, ethylene oxide sterilizing, or immersing the sensors in liquids.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

	Intravascular dyes such as indocyanine green or methylene blue may cause inaccurate readings.
	measurements may be adversely affected in the presence of higher ambient light.
	Excessive patient movement may cause inaccurate readings.
	The reusable patch is designed for 100 examinations. After this period, a possible degradation or malfunctions can occur.
	The device cannot be attached by disabled persons. Care-givers should perform the installation instead.
	If a serious adverse event occurs regarding our device, immediately report this situation to our company and to the competent authority of the European Union member state in which you reside.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	


Liability Notice

Failure to comply with the instructions outlined in this manual shall absolve Check Point R&D of any responsibility regarding the safety, reliability, and performance of this equipment. All operators must thoroughly read this manual before operating the system. Only authorized personnel may perform assembly, modification, or repairs. Electrical wiring must meet local regulatory standards. Use the equipment strictly in accordance with its intended use.

Terms of Warranty

Check Point R&D will, at its sole discretion, either repair or replace any part of its products found to be defective due to improper workmanship or materials. Repaired or replacement parts or products will be provided by Check Point R&D on an exchange basis. This warranty does not cover damage resulting from accidents, abuse or misuse, natural disasters or personal catastrophes, or unauthorized disassembly, modification, or repairs. CPC1S medical devices sold by Check Point R&D carry a warranty of 12 months from the date of purchase; all associated accessories, supplies, and disposable components carry a warranty period of 30 days from the date of purchase. This warranty is strictly limited to the repair or replacement of defective Check Point R&D products as described above. Check Point R&D shall not be liable for, and explicitly excludes from warranty coverage, any costs associated with patient care, servicing, or product installation. Check Point R&D commits to continued product support and will not discontinue or obsolete its products as long as necessary components and materials remain commercially available and reasonable customer demand persists.

Device Insulation Diagram



Manufacturer's EMC Declaration

Guidance and manufacturer's declaration – electromagnetic emissions		
This device is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an environment		
Emissions test	Compliance	Electromagnetic environment – guidance

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

RF emissions CISPR 11	Group 1	This device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

This device Image Intensifier is intended for use in the electromagnetic environment specified below. The customer or the user of This device Image Intensifier should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+ 6 kV contact (ECG Data) + 8 kV air (ECG Data)	+ 6 kV contact (ECG Data) + 8 kV air (ECG Data)	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz)magneticfield IEC 61000-4-8	3 A/m, 50Hz ECG	3 A/m, 50Hz ECG	If image distortion occurs, it may be necessary to position This device image intensifier further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

This device is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an electromagnetic environment.

IMMUNITY test	IEC 60601 TEST LEVEL	Compliance level	Electromagnetic environment – guidance
---------------	----------------------	------------------	--

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Conducted RF IEC 61000-4-6	150kHz-80MHz, 3V rms, 80% AM (1kHz) (6Vrms for ISM bands)	1 V _{active} 1 V _{active}	Portable and mobile RF communications equipment should be used no closer to any part of This device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter Recommended separation distance $d = 1.16\sqrt{P}$ $d = 1.2\sqrt{P} \text{ 0,15 MHz to 80 MHz}$ $d = 2.3\sqrt{P} \text{ 80 MHz to 2.7 GHz}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).^b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^c should be less than the compliance level in each frequency range.^d Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	80MHz - 2700MHz, 3V/m, 80% AM (1kHz)	10 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.



^b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.

^c Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the This device is used exceeds the applicable RF compliance level above, the This device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this device.

^d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 1 V/m.

Document No.	Release Date	Revision No.	Revision Date	USER MANUAL
TD.01.KK.02	15.03.2024	00	---	

Symbols and Meanings in the Label Prepared in Accordance with EN 15223-1 Standard

Guidance and manufacturer's declaration – electromagnetic immunity				
	Company Logo	Multiple Physiological Parameter Ambulatory Telemonitoring System Model CPC1S		Product Name
	Notified Body Number		Serial Number	
	UDI		Consult electronic instructions for use	
	Refer to instruction manual / booklet		Date of Manufacture	
	"CHECKPOINT R&D LTD" 37 Osvobojdenie Str, 6100 Kazanlak / Bulgaria		Keep Dry	
	Barcode Number		Type BF Applied Part	
	Medical Device		Class II equipment	
	Caution	IP22	IP Declaration	
	Latex Free		PVC Free	
	Bluetooth Wireless Communication Technology Not ionizing radiation – Device including Blue-tooth based RF transmitter.		Power Supply	
	≤93%, no condensation in storage		Storage Temperature: -20°C ~ 55°C	
TD.01.LT.01 / Rev.00 / ---				



37 Osvobojdenie str.
Kazanlak Bulgaria



info@checkpointcardio.com



0700 10 282